

# StratX™

Lung Analysis Platform

## TREAT WITH CONFIDENCE



Contact us for more information: [info@pulmonxstratx.com](mailto:info@pulmonxstratx.com)



[www.pulmonx.com](http://www.pulmonx.com)

**Pulmonx International Sàrl**  
Rue de la Treille 4  
2000 Neuchâtel, Switzerland  
Tel: +41 32 475 20 70  
Fax: +41 32 475 20 71

**Pulmonx Corporation**  
700 Chesapeake Drive  
Redwood City, CA 94063 USA  
Tel: +1 (650) 364-0400

#### REFERENCES:

<sup>1</sup> Klooster et al. N Engl J Med, 373;24, 2015. / <sup>2</sup> Koster TD, et al. Predicting Lung Volume Reduction after Endobronchial Valve Therapy Is Maximized Using A Combination of Diagnostic Tools. Respiration 2016 (in press). / <sup>3</sup> Reymond et al. AJR: 201, October 2013.

CAUTION: The Zephyr Endobronchial Valve is an investigational device in the US, limited by US law to investigational use.

CAUTION: The StratX Lung Analysis Platform is not available in the US.

©2016 Pulmonx Corporation. All Rights Reserved.

Pulmonx, Zephyr and StratX are registered or unregistered trademarks of Pulmonx Corporation.

# StratX™

## Lung Analysis Platform

A cloud-based quantitative CT analysis service that supports endobronchial valve (EBV) patient selection and treatment targeting by providing clinically-validated information on emphysema destruction, fissure completeness and volume

➤ Pulmonx EBV treatment is the most rigorously studied minimally-invasive treatment for severe emphysema and is proven to improve patients' breathing function, exercise capacity and quality of life.<sup>1</sup>

90%

### Accuracy

in identifying EBV responders and non-responders when combined with selective use of the Chartis Pulmonary Assessment System<sup>2</sup>




### Easy to Identify Patients Most Likely to Benefit

from EBV treatment



### Clinically-Validated and Consistent

across scanners and hospitals<sup>2</sup>



#### SUMMARY

Illustration summarizes key information

#### RESULTS

Table lists validated measurements by lobe:

- Fissure completeness
- Emphysema density (based on voxel density less than -910 HU)
- Inspiratory volume

# USER-FRIENDLY DESIGN

for clear interpretation and ease of use

# CLINICALLY-VALIDATED ALGORITHM

for fissure completeness





The StratX analysis quantifies the completeness of each fissure using an algorithm that has been validated in a retrospective study of over 200 EBV patients, the largest such analysis performed to date.<sup>2</sup> Fissure completeness is a proven predictor for volume reduction resulting from EBV therapy.<sup>3</sup>

RIGHT LUNG				
	RUL	RUL+RML	RML	RLL
% Fissure Completeness	64.2	65.0	84.9	65.0

## OPTIMAL APPROACH

The StratX analysis combined with selective use of the Chartis system results in higher accuracy than when either diagnostic tool is used alone.<sup>2</sup>

- 1 Assess all potential EBV patients with the StratX analysis.
- 2 For patients in which the StratX analysis identified partially complete fissures, proceed with a Chartis assessment.
- 3 Patients with complete fissures or identified with the Chartis assessment as CV- can receive EBV therapy.




StratX Analysis

**Complete Fissures**  
(>95% complete)

**Partially Complete Fissures**  
(80-95% complete)

**Incomplete Fissures**  
(<80% complete)



Chartis Assessment



EBV Therapy

CV-

CV+

**Do Not Treat with EBV Therapy**

# WORKFLOW WITH RAPID TURNAROUND TIME



## Capture CT Scan


Capture a high resolution chest CT scan according to the StratX CT parameters.



## Upload CT Scan

Use web browser to upload CT scan to the secure, cloud-based StratX platform.

- Automatically anonymized data with no patient health information transfer
- Secure 256 bit SSL socket level encryption



## Analyze Data + Generate Report

Data is analyzed by validated algorithms and the StratX report is uploaded to the StratX platform within 2-3 working days.



## Confidently Determine Treatment Options

Determine the most suitable treatment option for your patient using the quantitative StratX information and clinical judgment.



## Review Report

Access [pulmonxstratx.com](http://pulmonxstratx.com) to review the report in a .pdf (2D) or .html (3D) format from any clinical setting.